


Method for ECM Protein Extraction from Tissue



- 1) Pulverize tissue in liquid nitrogen
- 2) Homogenize tissue in detergent extraction buffer with BeadBeater or alternate homogenization tool.
- 3) Centrifuge sample at high speed, remove supernatant (Detergent Soluble Fraction) and proceed with protocol on remaining pellet.
- 4) Add additional detergent extraction buffer to pellet, vortex, spin, and combine supernatant with detergent soluble fraction from Step 3.
 - a. Repeat 1x.
- 5) Vortex remaining pellet with chaotrope extraction buffer, spin, and remove supernatant (Chaotrope Soluble Fraction). Proceed to step 6 with remaining pellet.
 - a. Repeat 1x.
- 6) Chemically digest remaining pellet in acidic environment for at least 18 hours.
- 7) Neutralize and wash pellet, solubilize in buffer, spin, remove supernatant (Chaotrope Insoluble Fraction).
 - a. Remaining pellet should contain less than 1% of total protein.
- 8) Proceed with digestion procedure -- We recommend FASP or GeLC-MS/MS for detergent/chaotrope removal.
 - a. Optional – Spike in Stable Isotope Label Peptides, QconCATs or AQUA, before or after digestion procedure, respectively, for downstream targeted MS analysis.
- 9) Extract Peptides and proceed with LC-MS/MS analysis for global protein identification or LC-SRM analysis for targeted protein ID and quantification.